




























structure	menref	name	loading	price	hazards
	H1396	Methylmorpholine on polystyrene, 0.8-1.0mmol/g	200-400 mesh 2% DVB	5g : £15.00, 25g : £60.00	Harmful
	H1411	Benzyl-diisopropylamine on polystyrene, 0.8-1.0mmol/g	200-400 mesh 2% DVB	5g : £15.00, 25g : £60.00	Air sensitive / Store under Nitrogen
	H1412	Methylthiomorpholine on polystyrene, 0.8-1.1mmol/g	200-400 mesh 2% DVB	1g : £8.00, 5g : £28.00, 25g : £110.00	Harmful
	H1413	Tris(2-aminoethyl)amine, polymer-supported, 0.6-1.0mmol/g	200-400 mesh 2% DVB	5g : £15.00, 25g : £60.00	Air sensitive / Store under Nitrogen
	H1415	4-Benzylpiperazine, polymer-supported, 0.8-1.1mmol/g on polystyrene	200-400 mesh 2% DVB	5g : £18.00, 25g : £70.00	Air sensitive / Store under Nitrogen
	H1416	Wang resin, 1% cross-linked, 0.8-1.0mmol/g [4-Benzylxybenzyl alcohol on polystyrene]	200-400 mesh	5g : £16.00, 25g : £65.00, 100g : £210.00	Harmful
	LH3	4-Benzylxybenzaldehyde polymer-supported, 0.8-1.0mmol/g on polystyrene	200-400 mesh 2% DVB	5g : £20.00, 25g : £75.00	Air sensitive / Store under Nitrogen
	LH4	S-Thiuronium chloride, polymer-supported, 0.8-1.0mmol/g on Merrifield resin [S-Methylisothiuronium chloride on polystyrene]	200-400 mesh 2% DVB	1g : £6.00, 5g : £22.00, 25g : £90.00	Harmful / Hygroscopic
	LH11	4-Nitrobenzoyl on polystyrene, 0.8-1.0 mmol/g	200-400 mesh 2% DVB	1g : £45.00, 5g : £180.00	Harmful
	LH13	p-Nitrophenylketoxime, polymer-supported, 0.8-1.0 mmol/g on polystyrene [Kaiser Oxime resin]	200-400 mesh 2% DVB	1g : £12.00, 5g : £42.00, 10g : £76.00	Harmful
	LH20	4-Hydroxythiophenol on Merrifield resin, 0.8-1.0 mmol/g	200-400 mesh 2% DVB	1g : £9.00, 5g : £30.00, 10g : £50.00	Irritant
	LH193	Wang resin, 2% cross-linked, 0.8-1.1mmol/g, 200-400 mesh [4-Benzylxybenzyl alcohol on polystyrene]	200-400 mesh 2% DVB	5g : £16.00, 25g : £65.00, 100g : £210.00	Harmful
	E3108	Benzaldehyde on polystyrene, 0.8-1.5 mmol/g [4-Formylpolystyrene]	200-400 mesh 2% DVB	1g : £12.00, 5g : £42.00, 10g : £76.00	Air sensitive
	E3109	Diphenylmethylphosphine, polymer-supported, 0.8-1.4mmol/g on polystyrene	100-200 mesh 2% DVB	1g : £6.00, 5g : £22.00, 25g : £90.00	Harmful

	E3110	Triphenylphosphine, polymer supported, 1.4-2.0mmol/g on polystyrene	200-400 mesh 2% DVB	5g : £24.00, 25g : £110.00, 100g : £400.00	Harmful
	LEH5	Bromopolystyrene, 1.0-2.0 mmol/g [Polystyrene, brominated]	200-400 mesh 2% DVB	5g : £10.00, 10g : £35.00	Harmful
	DH632	Benzenesulphonyl hydrazide, polymer-supported, 1.8-2.2 mmol/g on polystyrene [Sulphonyl hydrazide on polystyrene]	200-400 mesh 2% DVB	1g : £9.00, 5g : £30.00, 10g : £50.00	Harmful / Keep cold
	DH620	4-Benzyloxyacetophenone, polymer-supported, 0.8-1.1 mmol/g on Wang resin [Acetophenone on Wang resin]	200-400 mesh 2% DVB	1g : £12.00, 5g : £42.00, 10g : £76.00	Harmful
	DH621	4-Benzyloxybenzophenone, polymer-supported, 0.8-1.1 mmol/g on Wang resin [Benzophenone on Wang resin]	200-400 mesh 2% DVB	1g : £12.00, 5g : £42.00, 10g : £76.00	Harmful
	DH633	Benzeneboronic acid, polymer-supported, 2.6-3.2 mmol/g [4-Boronopolystyrene]	200-400 mesh 2% DVB	1g : £10.00, 5g : £40.00, 25g : £150.00	Harmful
	ER105	Benzylamine on polystyrene, 2.0-3.0 mmol/g [Aminomethylated polystyrene; AM resin]	200-400 mesh 2% DVB	5g : £30.00, 25g : £95.00	Air sensitive
	E3714	Methyl isothiocyanate, polymer-supported, 1.5-1.9mmol/g on polystyrene [Benzyl isothiocyanate on polystyrene]	200-400 mesh 2% DVB	1g : £12.00, 5g : £42.00, 10g : £76.00	Harmful / Moisture sensitive
	LH2	4-Benzyloxy-2,6-dimethoxybenzaldehyde, polymer-supported, 1.0-1.5mmol/g on polystyrene	200-400 mesh 2% DVB	1g : £12.00, 5g : £42.00, 10g : £76.00	Air sensitive
	ER441	Triphenylphosphine oxide, polymer supported, 1.2-1.8mmol/g on polystyrene	200-400 mesh 2% DVB	25g : £120.00, 100g : £400.00	Harmful
	GB15	N-Cyclohexylcarbodiimide, 0.8-1.0mmol/g on Merrifield resin	200-400 mesh 2% DVB	10g : £45.00, 25g : £90.00	Moisture sensitive / Store under Nitrogen
	E3894	Methyl isocyanate, polymer-supported, 1.5-1.9mmol/g on polystyrene [Benzyl isocyanate on polystyrene]	200-400 mesh 2% DVB	1g : £9.00, 5g : £30.00, 25g : £95.00	Harmful / Moisture sensitive
	LH15	p-Nitrophenol on benzylpolystyrene, 1.0-1.3mmol/g	200-400 mesh 2% DVB	1g : £20.00, 5g : £75.00, 10g : £120.00	Irritant